


## Workstation with thermal digestion unit



### Configuration

Thermal digestion unit (TAE) with controls integrated into the main operation unit

Alarm signal and automatic shutdown due to temperature monitoring

Ergonomically designed front

Convenient interchange from sliding doors to gloves or lifting gate

Horizontal ventilation for process optimization

Ventilation with microprocessor controlled system

System monitoring with sleep mode for downtime

H 14 Filter (HEPA, total retention > 99,995 %)

Are you interested in the wide range of combinations and possible solutions of this modular system?

Please contact us.

### Contact & Support

MK Versuchsanlagen und Laborbedarf e.K.

Stueckweg 10 - 35325 Muecke-Merlau

Tel.: +49 6400-9576030

Fax: +49 6400-9576031

info@mk-versuchsanlagen.de

www.mk-versuchsanlagen.de

