



- ▶ GITS®4 Software compliant with GAMP 5, CFR 21 Part 11 and EU Annex 11
- ▶ FDA and cGMP compliant
- ▶ ISO:9001 certified
- ▶ MK can help you to meet the Annex 1 glove leak testing requirements with minimum impact to your production schedules:

Integrity testing of the barrier systems, and leak testing of the glove system (...) performed minimum at the beginning and end of each batch is possible in an easy and fast way

To make your glove testing process easier and faster and for more information contact us:

**MK Versuchsanlagen und Laborbedarf e.K.**

Stueckweg 10  
35325 Muecke-Merlau  
Germany

sales@mk-versuchsanlagen.de  
☎ +49 6400 – 9576030  
www.mk-versuchsanlagen.de

*Scan me*  
↓



## Full Automatic Glove Tester

MK Glove Test Device  
Generation 3



Your solution to meet  
EU GMP Annex 1 changes!

- FDA and cGMP compliant
- Customizable to your port size
- 100% Made in Germany
- Easy cleaning, H<sub>2</sub>O<sub>2</sub>/VHP fumigable
- Complete qualification service
- Usable as stand-alone or fully integrated system




# Easy to Use

# Features

# Accessories and Services


## 1 Insert Glove Tester



### Total test time


- 8 – 15 Minutes (depending on glove material and port size)

Transportcase for smooth delivery (incl.)



### MK Glove Integrity Test System - GITS®

- Traceability of glove testing
- Port and glove recognition with integrated RFID technology
- Automatic machine validation and generated reports




Inductive battery (incl.) and charging station

### MK Glove Life Cycle System - GLCS®

- Monitoring of glove aging process
- Recognize glove failure in advance
- Glove life cycle tracking, from inventory to sterilization, through testing to replacement

Specific glove test parameters

Service trolley for up to 18 test devices with included charging stations





Qualification services (IQ, OQ, PQ and more)

Service and maintenance

WLAN site survey and planning

## 2 Turn on Glove Tester



## 3 Start measurement

